EnVision™ FLEX+, Mouse, High pH, (Link)

Code K8002

4th edition/ 4ème édition/ 4. Ausgabe

The kit contains reagents sufficient for 400-600 tests. For Autostainer Link instruments.

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Page/ Page/ Seite

ENGLISH

Intended Use	4
Summary and Explanation	4
Reagents	5
A. Materials provided	5
B. Materials required but not provided	5
C. Optional reagents	6
Precautions	7
Storage	7
Reagent Preparation	8
Specimen Collection and Preparation	8
Procedure	9
A. Pre-treatment procedure	9
B. Staining procedure	10
Quality Control	11
Interpretation of Results	11
Limitations	11
Explanation of symbols	30

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Intended Use

For in vitro diagnostic use.

Dako EnVision™ FLEX+ detection system is intended for use in immunohistochemistry together with Autostainer Link instruments. The system detects primary mouse and rabbit antibodies and the reaction is visualized by EnVision™ FLEX DAB+ Chromogen. If used with EnVision™ FLEX+ Mouse (LINKER) or EnVision™ FLEX+ Rabbit (LINKER) (Code K8009) signal amplification of primary mouse antibodies or primary rabbit antibodies, respectively, can be achieved. The EnVision™ FLEX+ reagents are intended for use on formalin-fixed, paraffin-embedded tissue sections.

Summary and Explanation

Dako EnVision™ FLEX+ detection system, has been designed to be a flexible system and gives an optimal staining on Autostainer Link instruments, when using the protocol recommended in this package insert. As a guideline, using EnVision™ FLEX+ Mouse (Code K8002) signal amplification increases **4-5 fold** and using EnVision™ FLEX+ Rabbit (Code K8009, optional reagent) increases signal amplification **2-3 fold**.

Prior to staining formalin-fixed, paraffin-embedded tissue sections should be subjected to deparaffinization and hydration followed by heat-induced epitope retrieval (HIER) using the target retrieval method specified in the package insert for the primary antibody. Please also see Procedure section for 3-in-1 specimen preparation. Code K8002 includes EnVision™ FLEX Target Retrieval Solution, High pH (50x). Alternatively, EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8004) or EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (Code K8005) is recommended. Some primary antibodies require enzymatic pre-treatment of tissue for optimal staining instead of HIER.

Endogenous peroxidase should be blocked with EnVision™ FLEX Peroxidase-Blocking Reagent (SM801) included in the kit. Due to an effective washing procedure and the presence of stabilizing proteins in the Dako reagents, extra blocking steps to reduce non-specific background staining are unnecessary.

EnVision™ FLEX Wash Buffer (20x) (DM831) included in the kit or available as optional reagent (Code K8007) is recommended.

Primary antibodies are not provided with the kit. We recommend the use of FLEX Ready-to-Use Primary Antibodies or Dako concentrated Primary Antibodies. EnVision™ FLEX Antibody Diluent (Code K8006) is recommended for the dilution of Dako concentrated Primary Antibodies.

EnVision™ FLEX+ Mouse (LINKER) (DM824) included in the kit or available as optional reagent (Code K8021) may be replaced by EnVision™ FLEX+ Rabbit (LINKER) (Code K8009) when using rabbit primary antibodies. EnVision™ FLEX+ Mouse (LINKER) and EnVision™ FLEX+, Rabbit (LINKER) may be applied for an optional signal amplification of mouse and rabbit primary antibodies respectively. As a guideline: EnVision™ FLEX+ Mouse (LINKER) provides 4-5 fold signal amplification and EnVision™ FLEX+ Rabbit (LINKER) provides 2-3 fold signal amplification.

Dako EnVision™ FLEX /HRP detection reagent (SM802) in the kit consists of a dextran backbone to which a large number of peroxidase (HRP) molecules and secondary antibody molecules have been coupled. A unique chemistry is used for the coupling reaction, which permits the binding of up to 100 HRP molecules and up to 20 antibody molecules per backbone.

The substrate system in the kit consists of two components: EnVision™ FLEX DAB+ Chromogen (DM827), a concentrated diaminobenzidine (DAB) solution, and EnVision™ FLEX Substrate Buffer (SM803) containing hydrogen peroxide. Before use EnVision™ FLEX DAB+ Chromogen must be diluted in EnVision™ FLEX Substrate Buffer. The substrate system produces a crisp brown end product at the site of the target antigen.

EnVision™ FLEX Hematoxylin (Code K8008) is recommended for counterstaining. The reagent provides a clear blue, nuclear staining.

The stained tissue sections may be mounted with either aqueous or organic-solvent-based mounting medium.

Reagents

A. Materials provided

SM801	EnVision™ FLEX PEROXIDASE-BLOCKING REAGENT EnVision™ FLEX Peroxidase-Blocking Reagent 3 x 40 mL, ready-to-use Phosphate buffer containing hydrogen peroxide, 15 mmol/L NaN ₃ and detergent.
SM802	EnVision™ FLEX /HRP EnVision™ FLEX /HRP 3 x 40 mL, ready-to-use Dextran coupled with peroxidase molecules and goat secondary antibody molecules against rabbit and mouse immunoglobulins. In buffered solution containing stabilizing protein and preservative.
DM827	EnVision™ FLEX DAB+ CHROMOGEN EnVision™ FLEX DAB+ Chromogen 3 x 3 mL 3,3'-diaminobenzidine tetrahydrochloride in organic solvent. The color of this reagent may vary from strong violet to colorless without having any influence on the performance of the kit.
SM803	EnVision™ FLEX SUBSTRATE BUFFER EnVision™ FLEX Substrate Buffer 12 x 20 mL Buffered solution containing hydrogen peroxide and preservative.
DM828	EnVision™ FLEX TARGET RETRIEVAL SOLUTION HIGH pH (50x) EnVision™ FLEX Target Retrieval Solution, High pH (50x) 9 x 30 mL, 50x concentrated Tris/EDTA buffer, pH 9
DM831	EnVision™ FLEX WASH BUFFER (20x) EnVision™ FLEX Wash Buffer (20x) 4 x 1 L, 20x concentrated Tris-buffered saline solution containing Tween 20, pH 7.6 (±0.1).
SM804	EnVision™ FLEX+ MOUSE (LINKER) EnVision™ FLEX+, Mouse (LINKER) 3 x 40 mL, ready-to-use Buffered solution containing stabilizing protein and an antimicrobial agent.

B. Materials required but not provided

Autostainer Link instrument

Dako PT Link

Dako Proteinase K, Ready-to-Use (if necessary)

Dako FLEX Ready-to-Use Primary Antibodies or suitably diluted concentrated primary rabbit or mouse antibodies from Dako

Dako FLEX Ready-to-Use Rabbit or Mouse Universal Negative Control or suitable negative control reagent for the primary antibody

Dako instrumentation utensils

Microscope slides, e.g. FLEX IHC Microscope Slides

General laboratory reagents for deparaffinization of paraffin-embedded tissue sections

Mounting medium (aqueous or organic-solvent-based) and coverslips

C. Optional reagents

K8004
DM828

EnVision™ FLEX
TARGET RETRIEVAL SOLUTION
HIGH pH (50x)

EnVision™ FLEX Target Retrieval Solution, High pH (50x)
3 x 30 mL, 50x concentrated
Tris/EDTA buffer, pH 9

K8005
DM829

EnVision™ FLEX
TARGET RETRIEVAL SOLUTION
LOW pH (50x)

EnVision™ FLEX Target Retrieval Solution, Low pH (50x)
3 x 30 mL, 50x concentrated
Citrate buffer, pH 6.1

K8006
DM830

EnVision™ FLEX
ANTIBODY DILUENT

EnVision™ FLEX Antibody Diluent
120 mL, ready-to-use
Tris buffer, pH 7.2, containing 15 mmol/L NaN₃, and protein.

K8007
DM831

EnVision™ FLEX
WASH BUFFER (20x)

EnVision™ FLEX Wash Buffer (20x)
1 L, 20x concentrated
Tris-buffered saline solution containing Tween 20, pH 7.6 (±0.1).

K8008
SM806

EnVision™ FLEX
HEMATOXYLIN

EnVision™ FLEX Hematoxylin
3 x 45 mL, ready-to-use
Aqueous solution of hematoxylin.

K8009
SM805

EnVision™ FLEX+
RABBIT (LINKER)

EnVision™ FLEX+ Rabbit (LINKER)
40 mL, ready-to-use
Buffered solution containing stabilizing protein and an antimicrobial agent.
Contains reagent sufficient for 130-200 tests based on the use of 200-300 µL per test.

K8021
SM804

EnVision™ FLEX+
MOUSE (LINKER)

EnVision™ FLEX+ Mouse (LINKER)
40 mL, ready-to-use
Buffered solution containing stabilizing protein and an antimicrobial agent.
Contains reagent sufficient for 130-200 tests based on the use of 200-300 µL per test.

Precautions

1. For professional users.
2. EnVision™ FLEX /HRP, EnVision™ FLEX+ Mouse (LINKER), and EnVision™ FLEX+ Rabbit (LINKER) contain material of animal origin and it cannot be excluded that trace amounts of human material could be present due to manufacturing procedures. As with any product derived from biological sources, proper handling should be used.
3. Do not expose EnVision™ FLEX Peroxidase-Blocking Reagent, EnVision™ FLEX /HRP, EnVision™ FLEX Substrate Buffer, EnVision™ FLEX DAB+ Chromogen, EnVision™ FLEX Substrate Working Solution or EnVision™ FLEX Hematoxylin to strong light during the procedure.
4. For EnVision™ FLEX Target Retrieval Solution, High pH (50x) and EnVision™ FLEX Target Retrieval Solution, Low pH (50x), a safety data sheet is available for professional users on request.
5. For EnVision™ FLEX Wash Buffer (20x), a safety data sheet is available for professional users on request.
6. For EnVision™ FLEX Peroxidase-Blocking Reagent, a safety data sheet is available for professional users on request.
7. EnVision™ FLEX DAB+ Chromogen contains 1–5% 3,3'-diaminobenzidine tetrahydrochloride, and is labeled:
Harmful.
R40 Limited evidence of a carcinogenic effect.
R43 May cause sensitisation by skin contact.
R68 Possible risk of irreversible effects.
S35 This material and its container must be disposed of in a safe way.
S36/37 Wear suitable protective clothing and gloves.
As a general rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper work procedure, the dangerous properties of the product and the necessary safety instructions (per European Union Directive 94/33/EC). Please refer to the Material Safety Sheet (MSDS) for additional information.
USA: 3,3'-diaminobenzidine (DAB) may be harmful by inhalation, in contact with skin and if swallowed. Material is irritating to eyes and skin. If skin contact should occur, rinse affected area with soap and water.
NOTE: Although diaminobenzidine is structurally related to benzidine there is no evidence for the carcinogenicity of diaminobenzidine. Consult Federal State, or local regulations for disposal.
8. EnVision™ FLEX Antibody Diluent (Code K8006) contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
9. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
10. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store EnVision™ FLEX and EnVision™ FLEX+ reagents at 2–8 °C.

Store EnVision™ FLEX Peroxidase-Blocking Reagent, EnVision™ FLEX /HRP and EnVision™ FLEX Substrate Buffer in the dark at 2–8 °C.

Store EnVision™ FLEX Hematoxylin (Code K8008) in the dark at room temperature. Do not use after expiration date stamped on the vials. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. The prepared Substrate Working Solution should be stored in the dark at 2–8 °C and used within 5 days. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact Dako Technical Services.

Reagent Preparation

A.1 EnVision™ FLEX Target Retrieval Solution, High pH (50x) (DM828, Code K8004)

Dilute a sufficient quantity of EnVision™ FLEX Target Retrieval Solution, High pH (50x) 1:50 using distilled or deionised water for the staining procedure that is planned. If Dako PT Link is used for pre-treatment, dilution can be performed by emptying the content of the Target Retrieval Solution (50x) vial into the Dako PT Link tank and adding distilled or deionised water to the marked line (Fill Line). Unused diluted solution may be stored at 2–8 °C for one month. Discard solution if cloudy in appearance. When used in PT Link for 3-in-1 specimen preparation procedure, the diluted solution can be used 3 times within a 5 day period if stored at room temperature.

A.2 EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (DM829, Code K8005)

Dilute a sufficient quantity of EnVision™ FLEX Target Retrieval Solution, Low pH (50x) 1:50 using distilled or deionised water for the staining procedure that is planned. If Dako PT Link is used for pre-treatment, dilution can be performed by emptying the content of the Target Retrieval Solution (50x) vial into the Dako PT Link tank and adding distilled or deionised water to the marked line (Fill Line). Unused diluted solution may be stored at 2–8 °C for one month. Discard solution if cloudy in appearance. When used in PT Link for 3-in-1 specimen preparation procedure, the diluted solution can be used 3 times within a 5 day period if stored at room temperature.

A.3 EnVision™ FLEX Wash Buffer (20x) (DM831, Code K8007)

Dilute a sufficient quantity of EnVision™ FLEX Wash Buffer (20x) 1:20 using distilled or deionised water for the staining procedure that is planned. Dilute by adding the concentrated Wash Buffer to a pre-measured quantity of distilled or deionised water to minimize foaming. Stir gently until the diluted solution appears homogenous. Unused diluted solution may be stored at 2–8 °C for one month. Discard solution if cloudy in appearance. When used in PT Link for 3-in-1 specimen preparation procedure, the diluted solution can be used 3 times within a 5 day period if stored at room temperature.

A.4 EnVision™ FLEX Substrate Working Solution

The DAB-containing EnVision™ FLEX Substrate Working Solution is prepared by mixing it thoroughly with 1 drop EnVision™ FLEX DAB+ Chromogen (DM827) per 1 mL EnVision™ FLEX Substrate Buffer (SM803). Use EnVision™ FLEX Substrate Working Solution within 5 days (store in the dark at 2–8 °C).

Specimen Collection and Preparation

The specimens may be formalin-fixed, paraffin-embedded tissue sections.

Fixation time is dependent on fixative and tissue type/thickness. For example, tissue blocks with a thickness of 3–4 mm should be fixed in neutral-buffered formalin for 18–24 hours.

The optimal thickness of paraffin-embedded sections is approximately 4 µm.

The specimens should be mounted on microscope slides, e.g. FLEX IHC Microscope Slides (Code K8020). The sections should be mounted on the slides as flat and wrinkle-free as possible. Too many wrinkles will have an impact on the staining results.

NOTE: The microscope slides must have a width suitable for the Autostainer Link instrument. Please refer to the Operator's Manual for the individual Dako instrument for definition of usable microscope slides.

Paraffin sections should be mounted from a pre-heated water bath containing distilled or deionized water. The water bath should contain no additives (such as gelatin, poly-L-lysine etc.). Sections should be dried by heating, generally at a temperature not above 60 °C for up to 60 minutes. To ensure proper adherence of sections to slides it is important to drain the water from beneath the sections prior to the oven drying process.

Procedure

Autostainer Link instruments use techniques based on different principles to obtain an optimal staining result. Before running protocols on your Autostainer Link instrument, please read carefully the **Operator's Manual for the dedicated Dako instrument**.

In order to ascertain reproducible epitope retrieval it is recommended always to load the slide holder fully with slides. This ensures an identical heating of sections in every run.

A. Pre-treatment procedure

Recommended 3-in-1 specimen preparation procedure using PT Link:

Deparaffinization, rehydration and heat-induced epitope retrieval (HIER) can be performed on formalin-fixed paraffin-embedded tissue sections using the 3-in-1 procedure:

1. Prepare a working solution by diluting the EnVision™ FLEX Target Retrieval Solution (50x) concentrate 1:50 in distilled or deionized water.
2. Fill PT Link tanks with sufficient quantity (1.5 L) of working solution to cover the tissue sections.
3. Set PT Link to pre-heat the solution to 65 °C.
4. Immerse the mounted, formalin-fixed, paraffin-embedded tissue sections into the pre-heated EnVision™ FLEX Target Retrieval Solution (working solution) in PT Link tanks and incubate for 20–40 minutes at 97 °C. The optimal incubation time should be determined by the user.
5. Leave the sections to cool in PT Link to 65 °C.
6. Remove each Autostainer slide rack with the slides from the PT Link tank and immediately dip slides into a jar/tank (e.g., PT Link Rinse Station, Code PT109) with diluted, room temperature EnVision™ FLEX Wash Buffer (20x).
7. Leave slides in the diluted, room temperature EnVision™ FLEX Wash Buffer (20x) for 1–5 minutes.
8. Place slides on an Autostainer Link instrument and proceed with staining. The sections should not dry out during the treatment or during the immunohistochemical staining procedure.
9. After staining it is recommended to perform dehydration, clearing and permanent mounting.

NOTE: When used in PT Link for 3-in-1 specimen preparation procedure, the diluted Target Retrieval Solution and diluted Wash Buffer can be used 3 times within a 5 day period if stored at room temperature.

HIER procedure using PT Link:

After deparaffinization and hydration to buffer (water), the tissue sections should be subjected to heat-induced epitope retrieval (HIER):

1. Prepare a working solution by diluting the EnVision™ FLEX Target Retrieval Solution (50x) concentrate 1:50 in distilled or deionized water.
2. Fill PT Link tanks with sufficient quantity (1.5 L) of working solution to cover the tissue sections.
3. Set PT Link to pre-heat the solution to 65 °C.
4. Manually deparaffinize and rehydrate tissue sections mounted on glass slides. For better adhesion of tissue to slide, use of coated or silanized slides is recommended.
5. Immerse room temperature tissue sections into the pre-heated EnVision™ FLEX Target Retrieval Solution (working solution), and incubate for 20–40 minutes at 97 °C. The optimal incubation time should be determined by the user.
6. Leave the sections to cool in PT Link to 65 °C.
7. Remove each slide rack from PT Link tank and immediately dip slides into a jar/tank (e.g., PT Link Rinse Station, Code PT109) with diluted, room temperature EnVision™ FLEX Wash Buffer (20x).
8. Leave slides in the diluted, room temperature EnVision™ FLEX Wash Buffer (20x) for 1–5 minutes.

9. Place slides on an Autostainer Link instrument and proceed with staining. The sections should not dry out during the treatment or during the immunohistochemical staining procedure.
10. After staining the slides should be mounted using aqueous or permanent mounting medium.

NOTE: When used in PT Link for HIER procedure, the diluted Target Retrieval Solution and diluted Wash Buffer can be used 3 times within a 5 day period if stored at room temperature.

HIER procedure using Coplin jars:

After deparaffinization and hydration to buffer (water), the tissue sections should be subjected to heat-induced epitope retrieval (HIER):

1. Place staining jars containing diluted EnVision™ FLEX Target Retrieval Solution (50x) in water bath.
2. Heat water bath and jars filled with EnVision™ FLEX Target Retrieval Solution (working solution) to 95–99 °C. Cover jars with lids to stabilize the temperature and avoid evaporation.
3. Immerse the room temperature sections in the preheated EnVision™ FLEX Target Retrieval Solution (working solution) in the staining jars.
4. Bring temperature of the water bath and EnVision™ FLEX Target Retrieval Solution back to 95–99 °C. Incubate for 20 (±1) minutes at 95–99 °C.
5. Remove the entire jar with slides from the water bath. Allow slides to cool in the EnVision™ FLEX Target Retrieval Solution for 20 (±1) minutes at room temperature.
6. Decant the EnVision™ FLEX Target Retrieval Solution and rinse sections with diluted, room temperature EnVision™ FLEX Wash Buffer (20x) for 1–5 minutes.
7. Place the rack with rinsed sections on an Autostainer Link instrument. Immediately begin the pre-programmed staining run.

NOTE: When used in Coplin jars HIER procedure, the diluted Target Retrieval Solution and diluted Wash Buffer can be used only 1 time (single use).

A few epitopes do not tolerate the HIER and some require enzymatic pre-treatment. Please refer to the package insert for the individual Dako primary antibody.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.

B. Staining procedure

Dako FLEX Ready-to-Use Primary Antibodies can be used with Dako EnVision™ FLEX+ detection system, applied on formalin-fixed, paraffin-embedded tissue sections.

A dilution guideline for Dako concentrated primary antibodies is provided in the package insert of the concentrated primary antibody.

The concentrated primary antibodies should be diluted in EnVision™ FLEX Antibody Diluent (DM830, Code K8006).

The staining steps and recommended incubation times are pre-programmed into the Autostainer Link software as the following visualization system protocols:

FLEX and FLEX 2x5 DAB (FLEX protocols)

FLEX+ Mouse and FLEX+ Mouse 2x5 DAB (FLEX+ Mouse (LINKER) protocols)

FLEX+ Rabbit and FLEX+ Rabbit 2x5 DAB (FLEX+ Rabbit (LINKER) protocols)

The recommended reagent application volume is 1 x 200 µL or 2 x 150 µL per slide. If the protocols are not available on the used Autostainer Link instrument, please contact Dako Technical Services.

Optimal incubation times of primary antibody and EnVision™ FLEX /HRP are dependent on the applied primary antibody. Please refer to the package insert for the individual Dako primary antibody. The user must verify the applied protocol.

When the staining procedure is completed, the specimens must be mounted. It is recommended to perform dehydration, clearing and permanent mounting when slides have been treated with the 3-in-1 specimen preparation procedure. Slides treated with HIER procedure can be mounted using aqueous

or permanent mounting medium. For aqueous mounting, mounting media such as Dako Glycergel™ Mounting Medium, Code C0563, or Faramount Aqueous Mounting Medium Code S3025 is recommended.

Quality Control

Each staining run should include a known positive control specimen to ascertain a proper performance of all the applied reagents. If the positive control specimen fails to demonstrate positive staining, labeling of test specimens should be considered invalid.

A negative control reagent should be used with each specimen to identify any non-specific staining. If non-specific staining cannot be clearly differentiated from the specific staining, the labeling of the test specimen should be considered invalid.

Interpretation of Results

The diaminobenzidine-containing Substrate Working Solution gives a brown color at the site of the target antigen recognized by the primary antibody. The brown color should be present on the positive control specimen at the expected localization of the target antigen. If non-specific staining is present, this will be recognized as a rather diffuse, brown staining on the slides treated with the negative control reagent. Nuclei will be stained blue by the hematoxylin counterstain.

Limitations

1. For automated use on Autostainer Link instruments only. When performing manual stainings, the optimal procedure must be determined and verified by the user, unless specified in the Package Insert for the specific Dako Primary Antibody
2. For use on formalin-fixed, paraffin-embedded tissue sections only. When using frozen sections or cell smears, the optimal procedure must be determined and verified by the user, unless specified in the Package Insert for the specific Dako Primary Antibody. Changes in the procedure may have impact on the number of tests provided with the kit.

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Dako Denmark A/S